

▶ Protein sizing by light scattering,  
molecular weight and polydispersity



Ulf Nobbmann <[Ulf.Nobbmann@malvern.co.uk](mailto:Ulf.Nobbmann@malvern.co.uk)>

Malvern Instruments

Malvern, Worcestershire WR14 1XZ

 [www.malvern.co.uk](http://www.malvern.co.uk)

© Malvern Instruments Ltd. - A **spectris** Company

# Outline

- ▶ Why light scattering?
  - Crystallization...
- ▶ Theory
  - Static light scattering (SLS) → molecular weight
  - Dynamic light scattering (DLS) → polydispersity
  - Electrophoretic light scattering (ELS) → zeta potential
- ▶ Application examples
  - Molecular weight
  - Sizing
  - Polydispersity
- ▶ Malvern Instruments & the Zetasizer Nano


## Why Light Scattering?

- ▶ The scattering intensity is a function of the molecular weight and concentration.
- ▶ Non-invasive technique, giving information on the size, mass, and charge of a protein sample.
- ▶ Light scattering is extremely sensitive to the presence of small amounts of aggregates.
- ▶ The velocity of a particle under an applied electric field is proportional to the charge.

# The direct link to Diamond

## ► Crystallization !

- Need crystal for structure determination
- Search for optimum conditions for growth
- Light scattering can give indication of likelihood of success





**Rule of thumb:**  
**Low polydispersity**  
**(Pd ≤ 20%) for best**  
**chance of crystal**  
**growth**

# Light - Matter Interactions


As light is sent through material there are several potential interactions:

- ▶ Transmission
- ▶ Absorption
- ▶ Fluorescence
  
- ▶ **Scattering !**



# Light - Matter Interactions: Scattering

The incident photon induces an oscillating dipole in the electron cloud. As the dipole changes, energy is radiated or scattered in all directions.




## Light Scattering

The scattering signal may be analysed by several methods:

- ▶ Average signal strength: **static**, 'classic'
- ▶ Fluctuations of signal: **dynamic**, quasi-elastic
- ▶ Shift of the signal: **electrophoretic**



# ▶ Static Light Scattering



## ▶ Molecular Weight Measurements




# Static Light Scattering (SLS)

Average scattering intensity is a function of the (particle) molecular weight and the 2<sup>nd</sup> virial coefficient.

## *Rayleigh Equation*

$$\frac{KC}{R_{\theta}} = \left( \frac{1}{M} + 2A_2C \right) \frac{1}{P(\theta)}$$



K = Optical constant  
M = Molecular weight  
 $A_2$  = 2<sup>nd</sup> Virial coefficient

C = Concentration  
 $R_{\theta}$  = Rayleigh ratio  
 $P(\theta)$  = Shape (or form) factor

# Static Light Scattering (SLS)

$$\frac{KC}{R_\theta} = \left( \frac{1}{M} + 2A_2C \right) \frac{1}{P_\theta}$$

$$K = \frac{2\pi^2}{\lambda_o^4 N_A} \left( n_o \frac{dn}{dc} \right)^2$$

$\lambda_o$  = laser wavelength  
 $N_A$  = Avogadro's number  
 $n_o$  = Solvent RI  
 $dn/dc$  = differential RI increment

$$P_\theta = 1 + \frac{16\pi^2 n_o^2 R_g^2}{3\lambda_o^2} \sin^2\left(\frac{\theta}{2}\right)$$

$R_g$  = Radius of gyration  
 $\theta$  = Measurement angle

$$R_\theta = \frac{I_A n_o^2}{I_T n_T^2} R_T$$


$I_A$  = Intensity of analyte (sample I – solvent I)  
 $n_o$  = Solvent RI  
 $I_T$  = Intensity of standard (toluene)  
 $n_T$  = Standard (toluene) RI  
 $R_T$  = Rayleigh ratio of standard (toluene)

## Static light scattering (SLS)

- ▶ The intensity of scattered light that a macromolecule produces is proportional to the product of the weight-average molecular weight and the concentration of the macromolecule ( $I \propto (M_w)(C)$ )
- ▶ For molecules which show no angular dependence in their scattering intensity, accurate molecular weight determinations can be made at a single angle (Rayleigh scatterers, isotropic scattering)
- ▶ This is called a Debye plot and allows for the determination of
  - ❑ Absolute Molecular Weight
  - ❑ 2nd Virial Coefficient ( $A_2$ )

## Debye plots: What do the measurements involve?

- ▶ Preparation of a number of concentrations of the unknown molecule (protein) in a suitable buffer



- ▶ Typical concentrations: 1, 2, 3 and 5 mg/mL

## Static Light Scattering (SLS)

$$\frac{KC}{R_{\theta}} = \left( \frac{1}{M} + 2A_2C \right) \frac{1}{P_{\theta}}$$

For Rayleigh scatterers,  $P(\theta) = 1$  and the equation is simplified to

$$\frac{KC}{R_{\theta}} = \left( \frac{1}{M} + 2A_2C \right) \quad (y = b + mx)$$

Therefore a plot of  $KC/R_{\theta}$  versus  $C$  should give a straight line whose intercept at zero concentration will be  $1/M$  and whose gradient will be  $A_2$

## Molecular Weight Example (Lysozyme in PBS)


$$\frac{dn}{dc} = 0.185(\text{mL} / \text{g})$$

$$I_{\text{tol}} = 192630 \text{ (counts/sec)}$$

$$I_{\text{sol}} = 21870 \text{ (counts/sec)}$$

Lysozyme Concentration (mg/mL)	Measured Intensity (counts/sec)	Intensity of Analyte (counts/sec)	KC/R <sub>0</sub> (1/Da)
1.006	87,830	65,960	6.1994 x 10 <sup>-5</sup>
3.018	222,900	201,030	6.4765 x 10 <sup>-5</sup>
5.029	366,770	344,900	6.6682 x 10 <sup>-5</sup>
10.059	742,570	720,700	6.7743 x 10 <sup>-5</sup>

# Molecular Weight Example (Lysozyme in PBS)



● KC/Ro                      — KC/Ro  
◇ Residual Intensity                      - - - Residual Intensity


## 2nd virial coefficient

- ▶ A thermodynamic property describing the interaction strength between the molecule and the solvent
- ▶ For samples where  $A_2 > 0$ , the molecules tend to stay in solution (protein molecules prefer contact with buffer)
- ▶ When  $A_2 = 0$ , the molecule-solvent interaction strength is equivalent to the molecule-molecule interaction strength – the solvent is described as being a theta solvent (protein doesn't mind buffer)
- ▶ When  $A_2 < 0$ , the molecule will tend to fall out of solution or aggregate (protein doesn't like buffer)





# ▶ Dynamic Light Scattering




## ▶ Molecular Size Measurements

# Dynamic Light Scattering

- ▶ Dynamic light scattering is a technique for measuring the size of molecules and nanoparticles
- ▶ DLS measures the time dependent fluctuations in the scattering intensity to determine the translational diffusion coefficient ( $D_T$ ), and subsequently the hydrodynamic radius ( $R_H$ )

# Dynamic Light Scattering

- ▶ Dynamic light scattering is a technique for measuring the size of molecules and nanoparticles
- ▶ DLS measures the time dependent fluctuations in the scattering intensity to determine the translational diffusion coefficient ( $D$ ), and subsequently the hydrodynamic size




- ▶ The rate of intensity fluctuation is dependent upon the size of the particle/molecule

# Dynamic Light Scattering (DLS)

Fluctuations are a result of Brownian motion and can be correlated with the particle diffusion coefficient and size.

*Stokes-Einstein*

$$R_H = \frac{kT}{6\pi\eta D}$$




q = Scattering vector  
 $R_H$  = Radius  
 T = Temperature

D = Diffusion coefficient  
 k = Boltzmann constant  
 $\eta$  = Solvent viscosity



# Brownian Motion

- ▶ **Random** movement of particles due to the bombardment by the solvent molecules that surround them




# Brownian Motion

- ▶ Temperature must be
  - accurately known (viscosity)
  - stable (otherwise convection present)
- ▶ The larger the particle the more slowly the Brownian motion will be
- ▶ The higher the temperature the more rapid the Brownian motion will be
- ▶ 'Velocity' of the Brownian motion is defined by the translational diffusion coefficient ( $D_T$ )

# Physical Constraints

The non-randomness of the intensity trace is a consequence of the physical confinement of the particles to be in locations very near to their initial locations across very short time intervals.




# Dynamic Light Scattering (DLS)

Fluctuations are a result of Brownian motion and can be correlated with the particle diffusion coefficient and size.

*Stokes-Einstein*

$$R_H = \frac{kT}{6\pi\eta D}$$



q = Scattering vector  
 $R_H$  = Radius  
 T = Temperature

D = Diffusion coefficient  
 k = Boltzmann constant  
 $\eta$  = Solvent viscosity



## Stokes-Einstein Equation

$$R_H = \frac{k_B T}{6 \pi \eta D_T}$$

where

$R_H$  = hydrodynamic diameter

$k_B$  = Boltzmann's constant

$T$  = absolute temperature


$\eta$  = viscosity

$D_T$  = diffusion coefficient


# Comparative Protein $R_H$ Values




**Lysozyme**  
 $M_W=14.5$  kDa  
 $R_H=1.9$  nm



**Insulin - pH 7**  
 $M_W=34.2$  kDa  
 $R_H=2.7$  nm




**Immunoglobulin G**  
 $M_W=160$  kDa  
 $R_H=7.1$  nm




**Thyroglobulin**  
 $M_W=650$  kDa  
 $R_H=10.1$  nm

5 nm


# Correlogram Interpretation




# Correlogram Interpretation




# Correlogram Interpretation



# Correlogram Interpretation




# Correlogram Interpretation



# Distributions By DLS

Comparison of Z average (Cumulants) size to multi-modal distribution results.



Z average: 12.4nm


(Single species assumption)



# Intensity, Volume And Number Distributions


Mixture containing equal numbers of 5 and 50nm spherical particles

NUMBER




VOLUME

$$= \frac{4}{3} \pi r^3$$



INTENSITY

$$= r^6$$




## Benefits Of Sizing By DLS

- ▶ Non-invasive
- ▶ High sensitivity ( $< 0.1$  mg/mL for typical proteins)
- ▶ Low volume (12  $\mu$ L)
- ▶ Scattering intensity is proportional to the square of the protein molecular weight, making the technique ideal for identifying the presence of trace amounts of aggregate.




# ▶ Electrophoretic Light Scattering



## ▶ Molecular Charge Measurements

# Electrophoretic Light Scattering (ELS)


Measured parameter is the frequency shift ( $\Delta\nu$ ) of the light scattered from a moving particle.



$\mu$  is the electrophoretic mobility,  $E$  is the electric field strength, and  $K$  is a constant.

# Measuring Electrophoretic Mobility

- ▶ Classical capillary electrophoresis (light microscope, stopwatch)



- ▶ The particles move with a characteristic velocity which is dependent on:


- Field strength
- Dielectric constant of medium
- Viscosity of the medium
- Zeta potential

# Electrophoresis

- ▶ Electrophoresis is the movement of a charged particle relative to the liquid it is suspended in under the influence of an applied electric field
- ▶ The electrophoretic mobility of a colloidal dispersion can be used to determine the zeta potential
- ▶ Zeta potential is the charge a particle acquires in a particular medium
- ▶ Zeta potential measurements can be used to predict dispersion stability
- ▶ Influenced by: pH, salts, concentration, additives,...

# Electroosmosis

## Electrophoresis in a Closed Capillary Cell



Electroosmosis is the movement of liquid relative to a stationary charged surface under the influence of an applied field

# Measuring Electrophoretic Mobility

- ▶ Laser Doppler electrophoresis (LDE)
  - Phase analysis light scattering (PALS)
  - Mixed mode measurements (M3)
- ▶ A laser beam is passed through the sample in the capillary cell undergoing electrophoresis
- ▶ Scattered light from moving particles is frequency shifted
- ▶ These small frequency shifts are measured
- ▶ The frequency shift  $\Delta f$  is equal to:

$$\Delta f = 2v \sin(\theta/2)/\lambda$$

$v$  = the particle velocity

$\lambda$  = laser wavelength

$\theta$  = scattering angle

...measure phase instead




## Mixed Mode Measurement (M3)


- ▶ Mixed mode measurement (M3) is a patented method that allows measurement at any point in a capillary cell
- ▶ It eliminates electroosmosis by reversing the applied field at a high frequency
- ▶ Malvern have combined M3 with PALS to improve the measurement sensitivity and accuracy (M3-PALS)

# Phase Analysis Light Scattering


## Phase Difference Demonstration




# Phase Plot From General Purpose



# Phase Plot From General Purpose




# Phase Plot From General Purpose



# Zetapotential Distribution Plot

## General Purpose




# Light Scattering Return

- Hydrodynamic Radius
- Distribution & Polydispersity
- Solution Composition
- Molecular Weight
- 2nd Virial Coefficient
- Conformation
- Shape Estimates
- Zeta Potential
- pI & Charge Estimates
- Formulation Stability




# ▶ Application Example



## ▶ Molecular Weight Measurements




# Absolute Protein Molecular Weight



# A<sub>2</sub> Crystallization Window

Crystallization Window:  $-0.8 > A_2 > -8 \times 10^{-4} \text{ mol mL} / \text{g}^2$ \*




*Rayleigh Equation*

$$\frac{KC}{R_{\theta}} = \frac{1}{M} + 2A_2C$$

\*George, A; Wilson, W.W. "Predicting protein crystallization from a dilute solution property", *Acta Crystallogr* **1994**, D50, 361-365.



# ▶ Application Example



- ▶
- ▶
- ▶
- ▶
- ▶

## Molecular Size Measurements


# Common Proteins

Comparison of known to estimated molecular weight for common proteins.

<b>Protein</b>	<b>Dia (nm)</b>	<b>Known MW (kDa)</b>	<b>Estimated MW (kDa)</b>	<b>%Error</b>
Lysozyme	3.8	14.7	15.1	-2.7
Chymotrypsinogen	4.8	25	26.1	-4.4
Carbonic Anhydrase	5.2	29	31.5	-8.6
Human Insulin (pH 7)	5.4	34.2	34.4	-0.6
Ovalbumin	6	43	44	-2.3
Hexokinase sub-unit	6.6	51	55	-7.8
Hemoglobin	7	65	63	3.1
Bovine Serum Albumin	7.1	67	65.3	2.5
Horse Alcohol Dehydrogenase	7.4	80	71.9	10.1
Amyloglucosidase	7.8	99	81.3	17.9
Hexokinase	8.6	102	102.2	-0.2
Yeast Alcohol Dehydrogenase	9.8	150	138.7	7.5
Apoferritin	16.4	443	462.7	-4.4
Thyroglobulin	20.2	669	753.5	-12.6

# Identifying Quaternary Structure

DLS results indicate an insulin structure that is dimeric at pH 2 and hexameric at pH 7, consistent with crystallographic data.



Human Insulin - pH 2  
 $R_H = 1.73$   
 Est MW = 12.1 kDa  
 Act MW = 11.4 kDa  
**Dimer**



Human Insulin - pH 7  
 $R_H = 2.69$   
 Est MW = 34.1 kDa  
 Act MW = 34.2 kDa  
**Hexamer**




# ▶ Application Example



## ▶▶▶▶ Polydispersity Measurements


# Polydispersity (Pd) From DLS

Pd is representative of the particle size distribution width.



*60 second measurement*

# Crystal Screening Using DLS






# Monoclonal Antibody Fragment


## Size, zeta potential and molecular weight

Size Distribution by Intensity




**Mean diameter = 5.1nm**

Zeta Potential Distribution



**Mean zeta potential = -7.6mV**



**$M_w = 20.7\text{KDa}$**   
 **$A_2 = - 0.0049 \text{ ml mol/g}^2$**




# ▶ Malvern Instruments




## ▶▶▶▶▶ & the Zetasizer Nano



# Zetasizer Nano



# Optics of the Zetasizer Nano



# Zetasizer Technical Specifications

Parameter	Value
Sizing range	0.6 nm to 6 $\mu$ m Diam
Concentration range	0.1 mg/mL Lys to 30w%
Min sizing sample volume	12 $\mu$ L
Min zeta sample volume	0.75 mL
Temperature control	2 to 90 $^{\circ}$ C
Conductivity range	0 to 200 mS/cm
Laser	3 mW 633 nm HeNe
Detector	APD



- ❖ Crystal screening
- ❖ Protein & polymer characterization
- ❖ CMC measurements
- ❖ Drug delivery systems
- ❖ Formulation stability
- ❖ Biological assemblies
- ❖ Virus & vaccine characterization
- ❖ Macromolecular critical points

## More information

- ▶ Application notes
- ▶ Multimedia presentations
- ▶ Brochures
- ▶ Detailed specifications

[www.malvern.co.uk/proteins](http://www.malvern.co.uk/proteins)